

Thermo Fisher Scientific Launches New Portable Electrochemistry Meter **Thermo Scientific Orion VERSA STAR Multi-Parameter Meter Platform Offers Up to Four Channels of pH/ORP, ISE, Conductivity and Dissolved Oxygen in Any Combination**

Dateline:

ORLANDO

Monday, March 12, 2012 9:17 am EDT

ORLANDO--([BUSINESS WIRE](#))--Thermo Fisher Scientific Inc., the world leader in serving science, today announced the new Thermo Scientific Orion [VERSA STAR](#) multi-channel, multi-parameter, bench-top meter platform. It offers up to four channels of pH/ORP, ISE, conductivity and dissolved oxygen in any combination. The company will showcase the new Thermo Scientific Orion VERSA STAR platform in Thermo Scientific booth 2665 during Pittcon 2012, being held from March 11-15, in Orlando.

“Modular design allows customers to get just the meter they need today while providing expansion capabilities to meet their needs tomorrow,” said Antonia Finlayson, product line manager, water analysis instruments, Thermo Fisher. “And a large, bright color display, multilingual interface, plain language prompts and custom display options make this meter extremely easy to use.”

The Orion VERSA STAR offers incremental methods for ISE, auto-recognition for polarographic or RDO optical dissolved oxygen probes and two Orion Star stirrer ports. Its 2000 memory datalog capacity, time/date-stamped datalogs and data transfer capability help labs meet GLP requirements. An RS232 communication port and two USB communication ports are standard.

The included electrode arm features a unique design that accommodates a variety of sensors at the exact height and angle needed to make an accurate measurement, even when using multiple probes of different types. Seven models of Orion [VERSA STAR](#) are available to meet single parameter or multi-parameter needs.

For more information about the new Thermo Scientific Orion [VERSA STAR](#), please visit Thermo Scientific booth 2665 at Pittcon 2012. Alternatively, please call 800-225-1480, email info.water@thermofisher.com or visit www.thermoscientific.com/OrionStarA.

For access to all Thermo Fisher news and product photos related to Pittcon 2012, please visit the online media room at www.thermofisher.com/news.

About Thermo Fisher Scientific

Thermo Fisher Scientific Inc. is the world leader in serving science. Our mission is to enable our customers to make the world healthier, cleaner and safer. With revenues of \$12 billion, we have approximately 39,000 employees and serve customers within pharmaceutical and biotech companies, hospitals and clinical diagnostic labs, universities, research institutions and government agencies, as well as in environmental and process control industries. We create value for our key stakeholders through three premier brands, Thermo Scientific, Fisher Scientific and Unity™ Lab Services, which offer a unique combination of innovative technologies, convenient purchasing options and a single solution for laboratory operations management. Our products and services help our customers solve complex analytical challenges, improve patient

“*Modular design allows customers to get just the meter they need today while providing expansion capabilities to meet their needs tomorrow*”

diagnostics and increase laboratory productivity. Visit <http://www.thermofisher.com>.

Photos/Multimedia Gallery Available: <http://www.businesswire.com/cgi-bin/mmg.cgi?eid=50200069&lang=en>

Contact:

Thermo Fisher Scientific Inc.
Sara Whyte, 978-232-6228
sara.whyte@thermofisher.com

<https://thermofisher.mediaroom.com/2012-03-12-Thermo-Fisher-Scientific-Launches-New-Portable-Electrochemistry-Meter>